


Hydrophobic interaction resins separates proteins based on differences in hydrophobicity, i.e., based on reversible interactions between proteins and hydrophobic groups on the surface of hydrophobic interaction resins. Hydrophobicity is enhanced at high ionic strengths and therefore binding in a high ionic strength environment is usually eluted by reducing the ionic strength. The unique adsorption separation mode makes hydrophobic interaction resins an ideal purification method after ammonium sulfate chromatography or after ion exchange high salt elution.




Product number	Product name	Spec	Ligand concentration $\mu\text{mol/mL}$	Particle size range μm	Maximum flow rate (cm/h)	Withstand pressure MPa	pH stability long-term [short-term]	Application characteristics	
HS060301025M		25mL							
HS060301100M		100mL							
HS060301500M		500mL							
HS060301001L	Phenyl Focurose FF (LS)	1L	20	45-165	400	≤ 0.3	3-13 [2-14]	Suitable for aromatic-containing aromatic ligands proteins	
HS060301005L		5L							
HS060301020L		20L							
HS060302025M		25mL							
HS060302100M		100mL							
HS060302500M		500mL							
HS060302001L	Phenyl Focurose FF (HS)	1L	40	45-165	400	≤ 0.3	3-13 [2-14]	High hydrophobicity and high loading capacity, suitable for aromatic ligands biomolecules	
HS060302005L		5L							
HS060302020L		20L							

Application Cases

Phenyl Focurose FF (HS) and Phenyl Focurose HP isolation of different hydrophobic proteins

Sample: 4mg/mL of mixed protein (Cytochrome C:
Ribonuclease A:Lysozyme=1:2:1)
Column: HT01, 1.0mL
Equilibrium solution: 0.1M Na₂HPO₄,
1.7M (NH₄)₂SO₄, pH7.0
Elution solution: 0.1M Na₂HPO₄, pH7.0
Flow rate: 1mL/min



≡ Tips

- ★ Different ligands and ligand concentration resins have different hydrophobic forces.
- ★ The salt concentration in the buffer is different for different protein hydrophobic interaction resins, or for purification using different hydrophobic interaction resins.
- ★ Temperature and pH have a great influence on protein hydrophobicity, and the pH and temperature should be constant during hydrophobic interaction resins.

Pre-assembled column ordering information

Product name	Spec	Product number	Product name	Spec	Product number
Phenyl Focurose FF (LS)	1mL	HS060301001E	Phenyl Focurose FF (HS)	1mL	HS060302001E
	5mL	HS060301005E		5mL	HS060302005E